# CompoSelect®

Blood component preparation with integrated leukocyte reduction



Blood component preparation with integrated leukocyte reduction



## Major benefits of CompoSelect® are:

- Full range of in-line filter systems suitable for different kinds of blood component preparation:
  - RCC in-line
  - RCC + Plasma in-line
  - RCC + Platelet in-line
  - WB in-line
- High leukocyte depletion acc. EU guidelines ( $<< 1 \times 10^6$  WBC/unit)
- Excellent physical characteristics for preparation and storage
- High quality, monitored assembly and welding for safe blood products
- User-friendly designed elements for efficient processing
- Well proven blood bag shape for high volume and cell recovery
- Wide compatibility with blood bag accessories and devices
- Fully implemented ISO norm 3826, part 1 3



### **Dimensions**

A variety of in-line blood bag systems with different sizes, additives and storage solutions are available for your individual requirements.

# **Ordering information**

PQ31450	CompoSelect® Quadruple CPD/SAG-M T&T WB, 450 ml, 15*/20** pcs/box
PQ31470	CompoSelect® Quadruple CPD/SAG-M T&T WB, 500 ml, 15*/20** pcs/box
PQ32250	CompoSelect® Quadruple CPD/SAG-M T&B RCC, 450 ml, 15*/24** pcs/box
PQ32270	CompoSelect® Quadruple CPD/SAG-M T&B RCC, 500 ml, 15*/24** pcs/box
PQ32350	CompoSelect® Quadruple CPD/SAG-M T&B RCC/Plasma, 450 ml, 15*/24** pcs/box
PQ31850	CompoSelect® Quadruple CPD/SAG-M - RCC/PLT, 450 ml, 15*/24** pcs/box

<sup>\*</sup> In Multivac packaging

These article codes are only exemplary. For more information, literature, technical details and working instructions as well as for equipment please contact your local sales representative or us.

#### CompoSelect® Premium

Р	S-single	1-CPDA-1	1-T&T	2-RCC	5 - 450 ml	0
Р	D-double	3-CPD/SAG-M	2-T&B	3-RCC/Plasma	7 - 500 ml	
Р	T-triple			4-WB		
Р	Q-quadruple					


Typical standard configurations are:

PQ31450 PQ32250 PQ32350 PQ31470 PQ32270 PQ31850 All articles with standard bag breaker

#### **Certifications**

CE marked Class IIb accord. EC-Directive 93/42/EEC

# Other products related to CompoSelect® that may be of interest to you



<sup>\*\*</sup> In Flow wrap packaging